Left ventricular restoration devices post myocardial infarction

Tom Hendriks¹ · Remco A. J. Schurer¹ · Lawien Al Ali¹ · Ad F. M. van den Heuvel¹ · Pim van der Harst¹

Published online: 17 May 2018 © The Author(s) 2018

Abstract

Even in the era of percutaneous reperfusion therapy, left ventricular (LV) remodeling after myocardial infarction (MI) leading to heart failure remains a major health concern. Contractile dysfunction of the infarcted myocardium results in an increased pressure load, leading to maladaptive reshaping of the LV. Several percutaneous transcatheter procedures have been developed to deliver devices that restore LV shape and function. The purposes of this review are to discuss the spectrum of transcatheter devices that are available or in development for attenuation of adverse LV remodeling and to critically examine the available evidence for improvement of functional status and cardiovascular outcomes.

Keywords Myocardial infarction · Left ventricular restoration · Transcatheter devices · Left ventricular remodeling

Introduction

Despite continuous improvements in acute and long-term treatment of myocardial infarction (MI) over recent decades, the resulting myocardial injury remains an important cause of adverse left ventricular (LV) remodeling, which over time progresses into the clinical syndrome of heart failure [1]. Typical for post-MI LV remodeling is dilatation, geometric change (from ellipsoid to spherical), and scar formation [2]. Both functional and structural measures of adverse LV remodeling, such as LV ejection fraction (LVEF) and LV dimensions, are associated with mortality and cardiovascular events after MI [3-6]. The beneficial effects of drugs or medical devices on LV remodeling measures have been associated with reduced long-term mortality [7]. In this review, we will focus on the role of percutaneous transcatheter delivered medical devices on LV remodeling. The definition of "medical device" varies according to local regulations, but generally, it is defined as an instrument, apparatus, software, or material that is intended for use in the diagnosis, treatment, monitoring, or prevention of disease [8]. It can be differentiated from a drug by its mechanism of action, which is not pharmacological, immunological, or metabolic. Numerous devices have been developed that attempt to restore the LV shape and revert adverse post-MI changes, intended to provide long-term benefits in terms of survival and functional status. Here, we will provide a critical overview of the empirical evidence supporting the use of LV restoration devices at various stages of LV remodeling post-MI, with a main focus on the transcatheter delivery route.

Remodeling after myocardial infarction

In the early stages after MI, influx of inflammatory cells leads to the destruction of collagen and regional thinning of the infarcted area [9]. The healing of necrotic tissue and formation of granulation tissue involves interactions among inflammatory cells such as macrophages, fibroblasts, and myofibroblasts [10]. As the infarcted region expands, cardiac fibroblasts deposit collagen in the infarct zone, which leads to scar formation and prevents further dilatation. Typical for post-MI remodeling is an increase in LV volumes and a more spherical LV shape due to a loss of myocytes and stretching of the remaining myocytes by the addition of sarcomeres in series [11]. Theoretically, dilatation can be beneficial by maintaining stroke volume through the Frank-Starling mechanism. However, the LV radius and the inverse of wall thickness are proportional with wall stress (law of Laplace). Consequently, LV dilatation increases wall stress and extends the burden on the remaining myocytes, leading to subendocardial myocardial ischemia and ultimately causing more damage [12]. Moreover, impaired LV contractility and reduced cardiac output after MI leads to the activation of neurohormonal



Pim van der Harst p.van.der.harst@umcg.nl

¹ University of Groningen, University Medical Center Groningen, Department of Cardiology, 9700RB Groningen, The Netherlands

pathways [13]. These pathways are thought to act as a way to maintain cardiac output through inotropic and chronotropic effects. However, these compensatory mechanisms result in an increased workload for the remaining myocytes, leading to progressive adverse remodeling. Drug therapies aimed at preventing heart failure post-MI mostly target neurohormonal pathways. Numerous studies have found evidence of a favorable effect of these drugs on parameters of LV remodeling and patient outcome [2]. Progressive dilatation is associated with larger infarct mass, nonscarred LV mass, and ongoing ischemia [14, 15]. It can lead to further hemodynamic consequences such as functional mitral regurgitation (FMR), which occurs in 20-25% of patients after MI and is associated with higher mortality rates [16]. Patients with progressive LV dilatation will eventually develop symptoms of heart failure. Despite very effective medical treatment options, the current chance of developing heart failure within 5 years after MI is 17% in men and 21% in women above the age of 45 [17], which occurs most distinctly in patients admitted with an anterior MI [18].

Guideline recommendations

Signs of LV remodeling can already be seen in the early stages after MI, within hours to days [2]. Contemporary guidelines recommend the use of transthoracic echocardiography (TTE) within the first 3 days after MI, to assess LVEF [19, 20]. A repeat TTE assessment 30-90 days after initial hospitalization is recommended because cardiac function can still recover in the case of myocardial stunning and hibernation. Cardiac magnetic resonance imaging (CMR) and computed tomography imaging (CT) are more suitable imaging modalities for research as a smaller sample size is needed to detect changes in LV structure, due to higher accuracy and reproducibility [21]. Recommended long-term therapies post-MI are based on cardiac rehabilitation, lifestyle interventions, and medical therapy. Contemporary MI guidelines do not recommend the use of transcatheter devices in long-term therapies post-MI [19, 20]. In the most recent heart failure guidelines [22–24], implantable cardioverter-defibrillators (ICDs), defibrillators with cardiac resynchronization therapy (CRT-Ds), and cardiac contractility modulation (CCM) are the only transcatheter devices recommended in the treatment of heart failure. This highlights that LV restoration devices have not yet been accepted as standard clinical practice.

Devices to reverse left ventricular remodeling

Short-term effects of drugs and devices on (reverse) LV remodeling are associated with reduced longer-term mortality rates [7]. Consequently, LV remodeling parameters are often used as alternatives or proxy variables for long-term mortality in clinical trials, because a smaller number of patients and a shorter follow-up are required to attain the same power. In the 1950s, the strong relation between adverse LV remodeling and patient outcome led to the hypothesis that surgical restoration of the original volume and ellipsoid shape of the LV could be beneficial in cases of severe LV remodeling. In the following paragraphs, we will discuss several surgical and transcatheter devices (Table 1) that have been investigated in post-MI patients in attempts to restore the LV shape as well as its hemodynamic and mechanical properties. The order of discussion is based on the indication for use.

Dilated ischemic cardiomyopathy

The first surgical attempt at LV restoration was the aneurysmectomy with a linear suture, first described by Cooley et al. in 1958 [25] and developed over the years [26]. The procedure was used in a selected group of patients with previous anterior MI and an aneurysmatic LV. The Dor procedure, using a circular suture and a pericardial patch that was covered by the residual myocardium, was later deemed superior [27]. In the STICH randomized controlled trial (N= 1000), CABG and surgical ventricular reconstruction reduced LV end-systolic volume by 19%, compared with 6% by CABG alone, but had no significant effect on mortality, hospitalization for cardiac events, and 6-min walk test (6MWT) distance [28]. A possible reason for this outcome proposed by the authors is that surgical reduction of LV volume, in addition to reducing wall stress, also reduces diastolic distensibility.

Several surgical LV restoration techniques have been investigated in patients with dilated cardiomyopathy, including not only ischemic but also idiopathic etiologies. In 1985, the first successful dynamic cardiomyoplasty procedure was performed, which was the surgical wrapping of an autologous latissimus dorsi muscle around the heart, which was activated by an external cardiomyostimulator for 10 weeks to gradually transform muscle fibers from type II to type I [29]. A prospective study including 68 patients observed a small increase in LVEF at 6 months (p = 0.05), but no significant change in peak oxygen consumption or cardiac index [30]. The cardiomyoplasty-skeletal muscle assist randomized trial (C-SMART) aimed to randomize 400 patients, but was terminated prematurely due to problems with patient recruitment and reimbursement. As opposed to autologous tissue, the LV can also be wrapped with synthetic material. Ventricular restraint therapy is the surgical placement of a multifiber polyester mesh around the LV, designed for patients with a dilated cardiomyopathy. It is intended to restore the ellipsoid shape of the LV and alleviate wall stress. In the ACORN trial (N = 300), the implanted device named the CorCapTM (formerly Acorn Cardiovascular, St. Paul, MN, USA) had no significant effect on 3- and 5-year mortality but did significantly reduce LV enddiastolic volume (LVEDV) up to 5 years after implantation, consistent across strata with and without mitral valve replacement [31–34]. It has to be noted that these results cannot be extrapolated to post-MI patients, since only 10% of the included patients had heart failure with an ischemic etiology.

Epicardial ventricular restoration

As an alternative to the invasive surgical ventricular reconstruction method, a minimally invasive surgical technique was designed to exclude the nonviable part of the LV in patients with LV dilatation and anteroseptal scarring. The ReviventTM myocardial anchoring system (BioVentrix, San Ramon, CA, USA), previously named the PliCath HFTM, is composed of polyester-covered titanium anchors (5 × 25 mm) mounted on a polyethylene-ether-ether-ketone tether, which are placed on the right side of the interventricular septum and on the LV wall. The anchors are drawn together to allow apposition of the LV free wall to the septum, thereby excluding the nonviable anteroseptal scar. Initial results in humans (*N*=11) demonstrate a stable reduction in LV volumes up to 12 months [35]. The Revivent system has received CE marking for commercial use in Europe.

Later developments of the device led to the introduction of a transcatheter component and a minor name change to Revivent-TC[™] (transcatheter) system, also receiving CE marking. The upgraded system is not completely transcatheter-based, still requiring a less invasive left thoracotomy. From outside the LV, a needle is used to puncture the LV wall, cross the LV, and puncture the interventricular septum. The needle position is monitored using fluoroscopy, and a Swan-Ganz catheter is introduced in the right jugular vein to monitor ventricular pressures. After reaching the right ventricle, the needle is replaced with a sheath and a guidewire. The guidewire is captured in the right ventricle by a snare from the Swan-Ganz catheter and connected with the internal anchor. The external anchor is placed from outside the LV wall. If necessary, additional pairs of anchors (2-3) are placed to achieve the desired line of apposition. Injection of contrast in the LV is used to confirm whether successful exclusion of the nonviable LV segment has been achieved. The Revivent-TC system has been successfully implanted in six sheep, successfully reducing LVESV, improving LVEF, and improving strain in border and infarct regions [36]. Study results from 51 patients treated with the Revivent and 20 patients treated with the Revivent-TC have been presented in the form of an abstract but have not yet been published in a peer-reviewed journal [37]. A prospective, multicenter, dual-arm pivotal study is currently aiming to include 146 subjects, who will be randomized to the Revivent-TC system or optimal medical therapy in a 2:1 ratio. Major inclusion criteria are the presence of an acontractile scar in the septal and anterior, apical or anterolateral regions of the LV, viable myocardium in the remote regions, LVEF under 45%, LVESVi larger than 50 ml/m², and NYHA class II or higher.

Transcatheter ventricular partitioning

The umbrella-like Parachute® device (Cardiokinetix, Redwood City, CA, USA) is intended to partition off the akinetic or aneurysmatic portion of the LV in patients with ischemic heart failure. The device is comprised of a fluoropolymer (ePTFE) membrane stretched over a selfexpanding nitinol frame, ranging between 65 and 95 mm in diameter when expanded. It is deployed into the LV apex and stabilized by 2-mm anchors at the end of each strut of the umbrella. The device provides efficacy by regional unloading of the akinetic LV region and global reduction of wall stress by reducing LV dimensions. Before implantation, LV anatomy has to be evaluated carefully, preferably by computed tomography imaging (CT), because anatomical characteristics such as prominent trabeculation or a "LV moderator band" are unsuitable for device implantation. The Parachute device has received CE marking for commercial use in Europe but is only approved for investigational use in the USA.

The device has been investigated in several PARACHUTE trials (Table 1). The PARACHUTE First-In-Human trial (N =39), which took place in Europe and the USA simultaneously, demonstrated safety and feasibility of the Parachute device in heart failure patients with LVEF between 15 and 40% and a dilated LV with an akinetic or dyskinetic anterior-apical wall [38]. There was a stable and significant reduction in LVEDV up to 3 years after implantation [39]. However, stroke volume and LVEF were also significantly lower at 3 years compared with measurements at baseline. The PARACHUTE III trial was a prospective, nonrandomized observational study conducted in Europe, including 100 subjects with similar inclusion criteria to the First-In-Human trial [40]. Procedural success was achieved in 97% of the subjects. LV volumes were significantly reduced at 12 months, and mean 6MWT distance improved from 372 m at baseline to 397 m (p < 0.01). During 12-month follow-up, cardiac mortality was 8.4%, devicerelated major adverse cardiac cerebral events (MACCE) occurred in 7.0%, heart failure hospitalizations in 24.1%, and stroke in 3% of subjects. Even though patients were given low-dose aspirin and warfarin for at least 12 months post device implant, the echo core lab identified thrombus on the device in 3.3% of the subjects at 12 months. The PARACHUTE IV was the first randomized controlled trial, aiming to include 478 patients with NYHA class III-IV, LVEF between 15 and 35%, LV motion abnormalities, and LV anatomy appropriate for Parachute implantation assessed by CT imaging and randomize to optimal medical therapy or Parachute implantation [41]. The primary endpoint was death or hospitalization for worsening heart failure. After including 331 subjects, the PARACHUTE IV trial was terminated in

June 2017 (clinicaltrials.gov \neq NCT01614652) and it is unclear whether investigation of the device will be continued.

In addition to the primary study results, there have been some interesting post hoc analyses. Hemodynamic assessment after Parachute implantation demonstrated an acute increase in stroke volume and mean aortic pressure [42]. A CT study investigating diastolic parameters before and after Parachute implantation observed that the Parachute device significantly reduced diastolic parameters such as the ratio between transmitral velocity in early diastole and early diastolic mitral septal tissue velocity (E/E'), and this reduction was correlated with improvement in NYHA functional class [43]. Another CT study observed a favorable effect of the Parachute on mitral valve geometry, by reducing interpapillary muscle distance and tenting height [44]. A pathology study investigated seven Parachute devices after cardiac transplantation (N=3)and at autopsy after deaths that were not device related (N = 4)[45]. The devices generally showed good healing by organized endocardial tissue growth on luminal and abluminal surfaces of the membrane. Parachute devices of long duration (>2 years, N=2) showed evidence of microscopic calcification, which could in theory lead to fatigue and rupture of the membrane on the long term.

The aforementioned study findings suggest that the Parachute device could be beneficial by reducing cardiac dimensions and end-diastolic wall stress and improving cardiac output. However, 3-year results demonstrate a reduction of LVEF and stroke volume. Results from the PARACHUTE IV trial will have to demonstrate the added benefit of the Parachute device in addition to optimal medical therapy. It is unclear whether the investigation and development of the device will be continued.

Ischemic mitral regurgitation

Several surgical LV techniques have been attempted to treat ischemic FMR before the use of a transcatheter delivery route. The Coapsys[™] (formerly Myocor, Maple Grove, MN, USA) device was designed to treat annular dilatation and papillary muscle displacement in patients with FMR by connecting an anterior and posterior epicardial pad with a subvalvular chord through the LV cavity and subsequently reducing the anteroposterior diameter by bearing load during the diastolic phase. The RESTOR-MV randomized controlled study was initiated, including patients referred for CABG with the presence of FMR [46]. Participants were stratified according to the need for mitral valve repair, as assessed by the surgeon. They were subsequently randomized to undergo either CABG (± mitral repair) or CABG and Coapsys implantation. The study was unfortunately terminated early due to insufficient funding. The investigators did find a significant survival advantage for patients who received the Coapsys device in the 165 subjects that were included, even though controls showed lower FMR grades during follow-up. In 2008, Edwards Lifesciences (Irvine, CA, USA) purchased the rights for the Coapsys device but since then has not revealed any plans for continuing the investigation or further development of the device.

Transventricular system

A transcatheter variant of the Coapsys, the iCoapsysTM repair system, was designed for the same patient category with ischemic FMR. The device consisted of the same elements as the original Coapsys (epicardial pads and a transventricular chord) but was implanted using transcatheter delivery. It was successfully implanted in 12 adult sheep, with no complications reported [47]. Prior to the start of a first-in-human study in 2008, Edwards Lifesciences purchased the rights for the iCoapsys device and discontinued the trial.

Similar to the iCoapsys repair system, the percutaneous septal sinus shortening system or PS^3 SystemTM has a bridge element inside the LV cavity, but instead of being attached on the epicardium, its anchors are placed on the right side of the atrial septum and in the great cardiac vein. It has been investigated in an ovine tachycardia model (N=19), successfully improving cardiac output and reducing FMR grade and septolateral diameter [48]. First-in-human implantations in two patients, immediately prior to explanation during planned clinically indicated surgical mitral valve repair, reduced septal-lateral dimensions and FMR grade [49].

Transvenous and subannular annuloplasty

Transcatheter interventions targeting the mitral apparatus such as percutaneous edge-to-edge mitral valve repair, percutaneous annuloplasty, and percutaneous mitral valve replacement are being reviewed elsewhere in this issue. A retrospective analysis of percutaneous edge-to-edge repair with the MitraClip (Abbott Vascular, Menlo Park, CA, USA) in 106 patients with ischemic FMR showed a significant reduction of left atrial volume at 1 year after implantation but no effect on LV volumes [50]. In the EVEREST II trial, MitraClip was associated with durable reduction of LVEDV up to 5 years, but this trial included mostly patients with degenerative mitral regurgitation [51]. In the following paragraph, we will discuss several transcatheter techniques here that target the coronary sinus system and the subannular space to reshape the LV with an effect on the mitral apparatus (Table 1).

A transcatheter mitral annuloplasty by implantation of a device in the coronary sinus (transvenous annuloplasty) has the advantage that no atrial septum puncture is required to gain access to the left heart. However, the coronary sinus is anatomically at a distance from the mitral annulus, and in some cases, the great cardiac vein passes over a coronary artery, which can be compressed or occluded by the device. The Monarc (previously Viking) system (Edwards Lifesciences, Irving, CA, USA) uses an implant with two self-expanding anchors and a spring-like "bridge," which was investigated in 72 patients of which 57% had a prior MI [52]. It was implanted in 59 patients (82%), and the primary safety endpoint of freedom from death, tamponade, or MI was 91% at 30 days and 82% at 12 months. A major issue in patients with implanted Monarc devices was coronary artery compression, which was observed in 15 patients who underwent followup angiography (30%), of which two presented with an acute MI. Another alternative, the percutaneous transvenous mitral annuloplasty (PTMA) device (formerly Viacor, Wilmington, MA, USA) consists of a PTMA catheter and nitinol rods of varying stiffness and is intended for patients with moderate to severe FMR, NYHA class II-IV, and LVEF 20-50%. The PTOLEMY-2 safety and feasibility trial included 43 subjects in whom implantation was attempted, resulting in 30 successful implantations of the device [53]. Unfortunately, four devices had to be removed at a later stage and two periprocedural deaths led to early termination of the trial by the sponsor. Also, during long-term follow-up after implantation, late erosions of the coronary sinus led to serious complications, in some cases fatal [54, 55]. The most thoroughly investigated transvenous annuloplasty device is the Carillon Mitral Contour System (Cardiac Dimensions, Sydney, Australia), which uses an implant that is composed of two self-expanding nitinol (nickeltitanium alloy) anchors connected by a curvilinear segment and delivered by a catheter via the right internal jugular vein, previously described in more detail [56, 57]. There have been three nonrandomized safety and efficacy studies conducted in Europe, including patients with at least moderate FMR, LVEF under 40%, NYHA functional class II or higher, and 6MWT distance between 150 and 450 m [56-58]. Where reported, FMR etiology was ischemic in 55-60% of patients. Study results demonstrated consistent low 30-day major adverse event rate (between 2 and 13%), significant reductions of regurgitant volume, and improvements in 6MWT distance. During follow-up evaluation in the TITAN trial, nine subjects (25% of subjects with a permanent implanted device) were observed to have a fractured anchor wire [57]. Although this was not linked to the occurrence of clinically relevant adverse events, the device (XE2) was modified to help reduce strain in the proximal anchor. The modified device (mXE2) was evaluated in the TITAN II trial, resulting in no more fractured devices [58]. The REDUCE FMR multicenter randomized controlled double-blind trial is currently randomizing 180 patients to the Carillon Mitral Contour System or optimal medical therapy in a 3:1 ratio, with primary efficacy endpoint being improvement in regurgitant volume at 12 months, assessed by an independent echocardiography core laboratory blinded to patient data [59]. A more recent device, the mitral loop cerclage (Tau-PNU Medical, Pusan, Korea), consists of a tension element made of stainless steel and an arch-like

coronary artery protection element. The device creates a loop through the coronary sinus and the right ventricle across the interventricular septum, covering the full circumference of the mitral annulus and applying circumferential tension. Because of the coronary artery protection device, it is suitable for a wider range of patients than the previously described devices. It has been investigated in a first-in-human feasibility study on five subjects with severe FMR and NYHA class III or IV, reporting a successful implantation in four subjects, reducing LVEDV and regurgitant volume [60]. Future studies will have to determine whether a transvenous mitral annuloplasty can be a feasible treatment for patients with ischemic FMR (Table 1).

The AccuCinch (Ancora Heart, Santa Clara, CA, USA) transcatheter direct mitral valve annuloplasty (or ventriculoplasty) system is designed to place anchors in the subannular space [68]. Cinching of these anchors reduces both basal LV and mitral annular dimensions. Anchors are delivered by anchor delivery catheters that are introduced through a modular guide tunnel via the femoral artery. A safety and feasibility trial is currently conducted, aiming to include a total of 40 patients.

Early phase after myocardial infarction

LV reconstruction techniques are intended to (partially) revert the changes to LV mechanical properties caused by myocardial injury. Ideally, this should be done during the early stages after MI, when scar formation, regional thinning of the myocardium, and LV dilatation can still be prevented. Most of the aforementioned studies, however, exclude patients within 90 days after MI. The reasoning behind this is that adverse LV remodeling can largely be reversed due to myocardial hibernation. Early assessment of myocardial viability by late enhancement CMR could be used to assess the risk of progressive LV remodeling and the need for an early intervention [69] (Fig. 1).

Biomaterial injection and infusion

As stated earlier, MI induces degradation of extracellular matrix, which reduces elasticity of the myocardium. A proposed option to restore the mechanical properties of the myocardium is by injecting biomaterials into the myocardium, which could serve as an artificial extracellular matrix. Biomaterial scaffolds incorporating stem cells were not included in this review, because we considered cellular products to be a different category. Several types of biomaterial injections have been tested in animal studies. The effect of calcium hydroxyapatite microsphere (CHAM) injections in infarcted myocardium was investigated in 31 sheep and appeared to limit LV dilatation and improve LVEF up to 4 weeks after ligation of the left anterior descending coronary artery [63]. A follow-up study

Interconserve to biotratical Sodium alginet (37%), strate alginet invessed metal algine vessed metal algine vessed metal vessed metal vessed metal vessed metal vessed metal vessed metal vessed metal vessed metal vessed metal vessed metal vessed metal vessed metal vessed metal vessed metal vessed ve		Type of device	Device	Indication*	Trials	Status	Primary endpoints	No. of subjects (nonischemic)
Fromatical guocome BCM reasonitarization RESERVATION I (e2) Completel Compl	arly phase	Intracoronary infusion of	Sodium alginate + calcium	2–5 days after large STEMI	First-in-human [61]	Completed	180-day symptomatic heart failure (3.7%), renal failure	27 (0)
Internocential Caloin hydroxygetis insection Decisity foot-Minute insection Animal statis (6:4-6) Completed Interception for the foot-Winute Minute Statis (7:4) 0 houtserial biochell instruct biochell instruct biochell instruct Biochell instruct Biochell instruct 0 houtserial biochell instruct biochell instruct Biochell instruct Biochell instruct 0 houtserial biochell instruct Biochell instruct Biochell instruct Biochell instruct 0 houtserial biochell instruct Biochell instruct Biochell instruct Biochell instruct 0 partitioning Partitioning Partitioning Completed Biochell instruct 0 partitioning Reviewer Display (2012) (201		biomaterial	gluconate BCM	and successful revascularization	PRESERVATION 1 [62]	Completed	(3.7%), stroke death (0%) 6-month ALVEDVi (<i>BCM 14.1</i> ± 28.9 ml/m ² vs. extinue 11 7 + 36.0 ml/m ² · n = 0.40)	306 (0)
nonuncut Injectible hydration Directly position Animal studies (6, 6) Completed		Intramyocardial injection	Calcium hydroxyapatite microspheres (CHAM)	Directly	Animal studies [63–65]	Completed	Infarct expansion, proteolytic pathways, $LVEF$, matrixed on the observation, MD conversion, MD convers	I
Moratial bydrogal (Ventricid) of ac STEM, bydrogal (Ventricid) Erst-in-human (XCT02350502) Rectining chonin SAE Ventricular partitioning Para-hut Completed Inplanation success and 6-month treedom fram device-claied MACE (745) Para-hut Completed Inplanation success and 6-month treedom fram device-claied MACE (745) Fortuning Revour-TC LVEF <456, ditatol I		of blomaterial	Injectable hyaluronic acid-based hydrogel	Directly post-MI	Animal studies [66, 67]	Completed	regional contractile strain, MK seventy Global LV remodeling, infarct thinning and expansion, and infarct stiffness	Ι
Ventricular Parachue Complete Implantation success and 6-month motion, suitable LV will motion, suitable LV will motion, suitable LV will motion, suitable LV antomy PARACHUTE Inst-in-human Completed Implantation success and 6-month metericular Frictioning Revivent-TC Revivent-TC Revivent-TC PARACHUTE V(CTD) 253335 Revining for versening heart matomy Frictioning Revivent-TC Revivent-TC Revivent-TC Revivent-TC Revivent-TC Frictioning Revivent-TC Revivent-TC Revivent-TC Revirention Revirention Frictioning Revivent-TC Revivent-TC Revirention Revirention Revirention Frictioning Revirenting Rev			Myocardial matrix hydrogel (VentriGel)	60 days–3 years after STEMI, 1 VFF 25–45%	First-in-human (NCT02305602)	Recruiting	6-month SAE	18 (0)
partonug motion, stable LV PARACHUTE III [40] Complete automy automonyconnous construction const	ate	Ņ	Parachute	$LVEF < 40\%$, dm_{model}	PARACHUTE fürst-in-human	Completed	Implantation success and 6-month	39 (0)
PARACHUTE IV [41] Terminated Death or rehospitalization for worsening heart inhue contractile scar. NYHA Revivent-TC LVEF < 45%, dilated LV, ReviventTC (<i>NCT0233378</i>) Recruiting Death or rehospitalization for worsening heart inhue contractile scar. NYHA Revivent-TC LVEF < 45%, dilated LV, ReviventTC (<i>NCT0233378</i>) Recruiting Chinis of the contractile scar. NYHA Revivent-TC LVEF < 45%, dilated LV, ReviventTC System	pitase			motion, suitable LV anatomy	PARACHUTE III [40]	Completed	1-year procedural- or device-related MACCE (77%)	100 (0)
Revivent-TC LVEF < 45%, dilated LV, acontractite sear, NYHA PARACHUTE V (NC70253535) Terminated could yo file Data contraction Iter and protection acontractite sear, NYHA Revivent-TC (NC70253535) Recuting class IL-IV Ecruing (Ericitronss, compared to surgical class IL-IV Compared to surgical Cospsys Severe FMR, NYHA class IL-III First-in-human (NC700512005) Terminated treatment and motical therapy class IL-IV Completed shortering (S) First-in-human (NC700512005) Precutaneous sepala sinus FMR Animal studies [48] Completed Sepal-lateral dimensions (- 17%), mean force on bidge (1, IC-I, 87 N) Monarc Moderate to severe FMR, mintal anuloplasty Erst-in-human (NC700512005) Terminated treatmet and motical therapy inter- and periprocedural Monarc Moderate to severe FMR, mintal anuloplasty Erst-in-human (NC700512005) Terminated treatmet and motical therapy inter of severe FMR, mintal anuloplasty ProLEMY-1/FE Nonarc Moderate to severe FMR, mintal anuloplasty Erst-in-human [99] Completed to the optical dimensions (- 29%), MR Or Spalateral dimensions (- 29%) MM 10,1% Spalateral dimensions (- 29%), MR Perturneous transvenous Moderate to severe FMR, mintal anuloplasty PrOLEMY-2[53]					PARACHUTE IV [41]	Terminated	Death or rehospitalization for worsening heart	331 (0)
Revivent-TC LVEF <45%, dilated LV, acontractile sear, NYHA Revivent TC (<i>NCT0253735</i>) Recuting Revivent FC System Factoring acontractile sear, NYHA Revivent TC System Coapsys Severe FMR, NYHA class II-III First-in-human (<i>NCT00512005</i>) Terminated Ferciveness, compared to surgical Coapsys Severe FMR, NYHA class II-III First-in-human (<i>NCT00512005</i>) Terminated Ferciveness, compared to surgical Percutaneous septal sinus FMR Animal studies [48] Completed Sepal-lateral dimensions (- 17%), mean shortening (P3 ³) First-in-human (<i>NCT00512005</i>) Terminated Sepal-lateral dimensions (- 17%), mean shortening (P3 ³) First-in-human (<i>NCT00512005</i>) Terminated Sepal-lateral dimensions (- 17%), mean shortening (P3 ³) First-in-human (<i>NCT0051205</i>) Terminated Sepal-lateral dimensions (- 17%), mean shortening (P3 ³) Moderate to severe FMR, EVOLUTION [52] Completed Sepal-lateral dimensions (- 17%), mean miral amuloplasty NYHA class II-IV, LVEF EVOLUTION [53] Completed Sepal-lateral dimensions (- 17%), mean filated LV, LVEF State Terminated Sepal-lateral dimensions (- 17%), mean Sepalateral dimensions (- 17%), mean					PARACHUTE V (<i>NCT02543632</i>)	Terminated	tatiure Onality of life	85 (0)
Image: Interpret and the servent TC System Recruiting contractile sear, NYHA Revirent-TC System Recruiting consest compared to surgical vector of class II-IN class II-III First-in-human (<i>NCT00512005</i>) Terminated Effectiveness, compared to surgical vector of class II-IN class II-III First-in-human (<i>NCT00512005</i>) Terminated Restinent and medical therapy (<i>NCT0291240</i>) Percutaneous septal sinus FMR Aimal studies [48] Completed Severe FMR, <i>NFHA</i> class II-III First-in-human (<i>NCT00512005</i>) Terminated Restinent and medical therapy (<i>NCT02912405</i>) Percutaneous septal sinus FMR Aimal studies [48] Completed Septal-lateral dimensions (- <i>17%</i>), mean (<i>not completed dilated LV, LVEF</i> Septal-lateral dimensions (- <i>106-N</i>), mean (<i>not completed dilated LV, LVEF</i> Monarc Moderate to severe FMR, intra anuloplasty NYHA class II-IV, LVEF Completed 30-day freedom from dath, tamponde, and MI (<i>97%, 82% at 1 year</i>) Percutaneous tansverous Moderate to severe FMR, intra anuloplasty PTOLEMY-2 [53] Completed 30-day freedom from MACE3 and dilated LV, LVEF Percutaneous tansverous Moderate to severe FMR, intra anuloplasty MODEMY-2 [55] Completed 30-day freedom from MACE3 and dilated LV, LVEF Percutaneous tansverous Moderate to severe FMR, intra NFS Terminated 30-day freedom from MACE3 and dilated LV. LVEF 30-day freedom from MACE3 and dilated LV. LVEF Percutaneous tansverous Moderate to severe FMR, intrean		Epicardial	Revivent-TC	LVEF < 45%, dilated LV,	ReviventTC1 (NCT02553785)	Recruiting	6-month SAE and ALVESVi	50(0)
Coapesys class II-IV class II-II Connected and operational study (NCT00512005) Terminated and moderation sefery and MR reduction and MR reduction Percutaneous septal sinus FMR Animal studies [48] Completed Septal-lateral dimensions (- 17%), mean affery and MR reduction Percutaneous septal sinus FMR Animal studies [48] Completed Septal-lateral dimensions (- 17%), mean affery and MR reduction shortening (PS ¹) First-in-human (NCT00512005) Terminated Septal-lateral dimensions (- 17%), mean affery and MR reduction Monarc Moderate to severe FMR, dilated LV, LVEF EVOLUTION [52] Completed 30-day freedom from death, tamponade, and MI (9/%, 82% at 1) sev) Perculaneous transvenous Moderate to severe FMR, mirral annuoplasty PTOLEMY-2 [53] Terminated Perculaneous transvenous Moderate to severe FMR, mirral annuoplasty AMDEUS [56] Completed 30-day freedom from MACEs and 6-month MR reduction (PTMA) 20-50%, mirral annuoplasty Trannited 30-day MAE (13%) 6-month MR reduction (PTMA) 20-50%, mirral annuoplasty Trannited 30-day MAE (13%) 6-month MR reduction (PTMA) 20-30%, mirral annuoplasty MADEUS [56] Completed 30-day MAE (13%) 177N1/57 (PTMA) 20-30%, mirral annuoplasty MADEUS [59] Completed 30-day MAE		ventricular restoration	_	acontractile scar, NYHA	Revivent-TC System	Recruiting	Effectiveness, compared to surgical	126 (0)
Percutaneous septal sinus FMR Animal studies [48] completed safety and MR. reduction shortening (P3 ³) First-in-human [49] First-in-human [49] Completed Septal-lateral dimensions (- 17%), mean force on bridge (1.16-1.87 M) system Monarc Moderate to severe FMR, initial annuloplasty EVOLUTION [52] Completed Septal-lateral dimensions (- 187, M) sty Monarc Moderate to severe FMR, initial annuloplasty EVOLUTION [52] Completed 30-day freedom from death, tamponade, and M1 0/6, 82% at 1 year) sty Percutaneous transvenous Moderate to severe FMR, initial annuloplasty MADEUS [56] Completed 30-day freedom from MACEs and 6-month MR reduction of Ordino from MACEs MADEUS [56] Completed 30-day MAE (13%) of Ordino from MACEs and 6-month MR reduction 0-day MAE (13%) of Moderate to severe FMR, intral loop cerclage MAC (13%) Mitral loop cerclage Mitral loop cerclage First-in-human [60]		Transventricular system	iCoapsys	Severe FMR, NYHA class II–III	First-in-human (NCT00512005)	Terminated	ucauncil and incurcal uncapy Intra- and periprocedural	j (?)
Minimation Animation Animation Animation system system space sequence Space and state dimensions (-1, 1/6-1/37, N), mean system system monarc Moderate to severe FMR, dilated LV, LVEF > 25% EVOLUTION [52] Completed Septal-lateral dimensions (-29%), MR grade reduction (1 <i>and</i> 2) System monarc Moderate to severe FMR, dilated LV, LVEF EVOLUTION [52] Completed Septal-lateral dimensions (-29%), MR grade reduction (1 <i>and</i> 2) System mithal ammloplasty NYHA class II-V, LVEF PTOLEMY-2 [53] Terminated 30-day freedom from MACEs and 6-month MR reduction Ordsy freedom from MACEs and (PTMA) 20-50% AMADEUS [56] Completed 30-day MAE (13%) Ordsy freedom from MACEs and 6-month MR reduction 0-day MAE (13%) 6-month MR reduction (PTMA) 20-50% TTAN II [58] Completed 30-day MAE (13%) (PTMA) 20-50% TTAN II [58] Completed 30-day MAE (13%) (PTMA) 20-50% TTAN II [58] Completed 30-day MAE (13%) (PTMA) 20-50% TTAN II [59] Completed 30-day MAE (13%) (PTMA) 20-50% TTAN II [59] Completed 30-day MAE (21%) (PTMA) 6/month MR reduction 1-yeerdon from MACE 9/mono			Domitor for the state of the st		المسابعة المسابعة	Commission	safety and MR reduction	
system system system first-in-human [49] Completed Septal-lateral dimensions (- 29%), MR grade reduction (1 and 2) sty Monarc Moderate to severe FMR, dilated LV, LVEF > 25% EVOLUTION [52] Completed Septal-lateral dimensions (- 29%), MR grade reduction (1 and 2) sty Percutaneous transvenous Moderate to severe FMR, multiplasty NVHA class II-IV, LVEF PTOLEMY-2 [53] Terminated 30-day freedom from death, tamponade, and MI (97%, 82% at 1 year) opTMA) 20-50% TTAN [57] Completed 30-day freedom from MACEs and 6-month MR reduction nMMA) 20-50% TTAN [57] Completed 30-day MAE (13%) off Moderate to severe FMR, moderate to severe for the			rercutaneous septat sinus shortening (PS ³)	FINIK	Animai studies [48]	Completed	Septal-lateral undensions (-1.7%) , mean force on bridge $(1.16-1.87 N)$	I
Monarc Moderate to severe FMR, dilated LV, LVEF > 25% EVOLUTION [52] Completed 30-day freedom from death, tamponade, and MI (91%, 82% at 1 year) Percutaneous transvenous Moderate to severe FMR, mitral amuloplasty PTOLEMY-2 [53] Terminated 30-day freedom from MACEs and 6-month MR reduction PTMA) 20–50% TTAN [57] Completed 30-day freedom from MACEs and 6-month MR reduction Moderate to severe FMR, (PTMA) Moderate to severe FMR, aliated LV, LVEF < 50%, dilated LV, LVEF < 50%, mitran list AMADEUS [56] Completed 30-day MAE (1.9%) Moderate to severe FMR, for the from MACEs and dilated LV, LVEF < 50%, dilated LV, LVEF < 50%, mitran list AMADEUS [56] Completed 30-day MAE (1.9%) Moderate to severe FMR, for the from death, tamponade, dilated LV, LVEF < 50%, for the from MACEs and for the from MAE (1.9%) So day MAE (1.9%) So day MAE (1.9%) Mitral loop cerclage Mitral loop cerclage Terminated 30-day MAE (1.9%) So day MAE (1.9%) Mitral loop cerclage Mitral loop cerclage Terminated 1-year freedom from MAE. A regurgitant volume, and clinical composite of death, time to first HF hospitalization			system		First-in-human [49]	Completed	Septal-lateral dimensions (-29%) , MR	2 (1)
Percutancous transvenous mitral amuloplastyPercutancous transvenous Moderate to severe FMR, $20-50\%$ PrOLEMY-2 [53]Terminated 0 -month MR reduction 0 -month MR reductionPTIMA) $20-50\%$ 0 -month MR reduction 0 -day freedon from MACEs and 0 -month MR reductionPTIMA) $20-50\%$ 0 -month MR reduction 0 -day MAE (13%) 0 -day MAE (13%)CarillonModerate to severe FMR, dilated LY, LVEF < 50% , 0 -MWT distance $150-450$ m REDUCE FMR [59]Completed $30-day MAE (1.9\%)0-day MAE (1.9\%)More that the severe FMR,dilated LY, LVEF < 50\%MAE (1.9\%)0-day MAE (1.9\%)0-day MAE (1.9\%)Moderate to severe FMR,dilated LY, LVEF < 50\%,0-MWT distance 150-450 mREDUCE FMR [59]Completed30-day MAE (1.9\%)0-day MAE (2.8\%)More that the severe transmitfor the teach from the severe transmitreconting1-year A regurgitant volume,recontingMitral loop cerclagefirst-in-human [60]CompletedMitral loop cerclagefirst-in-human [60]Completed$		Transvenous	Monarc	Moderate to severe FMR, dilated TV TVFF ~ 35%	EVOLUTION [52]	Completed	grade reduction (1 and 2) 30-day freedom from death, tamponade, and MI (01%, 82%, at 1 your)	72 (23)
Moderate to severe FMR, dilated LV, LVEF < 50%, MWT distance 150-450 m TITAN [57] Completed 30-day MAE (13%) 6MWT distance 150-450 m TITAN II [58] Completed 30-day MAE (2.8%) 6MWT distance 150-450 m TITAN II [58] Completed 30-day MAE (2.8%) 6MWT distance 150-450 m TITAN II [58] Completed 30-day MAE (2.8%) 6MWT distance 150-450 m TITAN II [58] Completed 30-day MAE (2.8%) 6MWT distance 150-450 m REDUCE FMR [59] Active, not 1-year A regurgitant volume recontiting 1-year freedom from MAE, A regurgitant volume 1-year freedom from MAE, A regurgitant volume control 1-year freedom from MAE, A regurgitant volume 1-year freedom from MAE, A regurgitant volume control 1-year freedom from from MAE, A regurgitant volume 1-year freedom from MAE, A regurgitant volume control 1-year freedom from from MAE, A regurgitant volume 1-year freedom from MAE, A regurgitant volume control 1-year freedom from from from from from from from fr		Constantin	Percutaneous transvenous mitral annuloplasty (PTMA)	Moderate to severe FMR, NYHA class II–IV, LVEF 20–50%	PTOLEMY-2 [53]	Terminated	30-day freedom from MACEs and 6-month MR reduction	43 (13–26)
6MWT distance 150–450 m TITAN II [38] Completed 30-day MAE (2.8%) REDUCE FMR [39] Active, not 1-year Δ regurgitant volume certaiting 1-year freedom from MAE, Δ regurgitant CARILLON (<i>NCT03142152</i>) Recruiting 1-year freedom from MAE, Δ regurgitant volume, and clinical composite of death, time to first HF hospitalization and improvement in 6MWT distance First-in-human [60] Completed			Carillon	Moderate to severe FMR, dilated LV, LVEF < 50%,	AMADEUS [56] TITAN [57]	Completed Completed	30-day MAE (13%) 30-day MAE (1.9%)	48 (?) 53 (19)
CARILLON (NCT03142152) Recruiting 1-year freedom from MAE, A regurgitant volume, and clinical composite of death, time to first HF hospitalization and improvement in 6MWT distance First-in-human [60] Completed				6MWT distance 150-450 m	TITAN II [58] REDUCE FMR [59]	Completed Active, not	30-day MAE (2.8%) 1-year Δ regurgitant volume	43 (15) 180 (?)
First-in-human [60] Completed					CARILLON (NCT03142152)	Recruiting	1-year freedom from MAF, A regurgitant volume, and clinical composite of death, time to first HF hospitalization	400 (?)
			Mitral loop cerclage		First-in-human [60]	Completed	and improvement in olviw 1 distance	5 (2-4)

 $\underline{\textcircled{O}}$ Springer


Type of device	Device	Indication*	Trials	Status	Primary endpoints	No. of subjects (nonischemic)
		Severe FMR, NYHA class III-IV			1-month freedom from MACE (80%), reduction	
Direct mitral valve annuloplasty (ventriculoplasty)	AccuCinch	Moderate to severe FMR, NYHA class II-IV, LVEF ≥20%	Safety and efficacy trials (<i>NCT00800046</i> , <i>NCT01899573</i> , <i>NCT02806570</i> ,	Recruiting 70,	ou regurgitant volume (00.5%) and encouve regurgitant orifice area (76%) 30-day device-related or procedure- related MAE, MACE, acute and 30-day MR reduction	197 (?)
			NCT02153892, NCT03183895)			

diastolic volume indexed to body surface area; LVESV, left ventricular end-systolic volume indexed to BSA; MACCE, major adverse cardiac cerebral events; MACE, major adverse cardiac events; MAE, major adverse events; MI, myocardial infarction; MR, mitral regurgitation; No., number; NYHA, New York Heart Association; SAE, serious adverse events; STEMI, ST-elevation myocardial infarction *Inclusion criteria occasionally vary in case of multiple trials. The indications mentioned below are generalized and not intended to reflect the full scope of in- and exclusion criteria randomized 24 adult male sheep to receive either saline or CHAM injection after ligation of the left anterior descending coronary artery [64]. In addition to improving LVEF and reducing LVEDV compared to controls, the authors observed a reduction of collagen content in the sheep that received CHAM injections. A more recent study randomized 25 Yorkshire swine to receive either 20 normal saline or CHAM injections at 30 min after ligation of the circumflex artery [65]. The animals receiving CHAM injections showed reduced infarct thinning and progressive improvement in border zone and global LV function. To date, no studies in humans have been reported.

Another widely investigated biomaterial is injectable myocardial matrix hydrogel, which is produced by decellularizing porcine myocardial tissue and processing it to form a myocardial matrix with the ability to gel upon injection. Myocardial matrix hydrogel injections have been investigated in 58 rats undergoing ischemia reperfusion followed by injection of either the hydrogel or saline 2 weeks later [70]. The hydrogel increased endogenous cardiomyocytes in the infarct area and maintained cardiac function without inducing arrhythmias. Similar results were observed in a study on 10 pigs, using transcatheter injections [71]. The authors also investigated the hemocompatibility of porcine myocardial matrix with human blood and observed no effect on clotting times. Possible mechanisms underlying the benefit of myocardial matrix injections are an altered inflammatory response, reduced cardiomyocyte apoptosis, enhanced neovascularization, diminished cardiac hypertrophy and fibrosis, and enhanced recruitment of progenitor cells [72]. A phase 1 study is currently enrolling post-MI patients to study the safety of transcatheter administration of myocardial matrix hydrogel in humans.

A less frequently investigated type of hydrogel that has been investigated in animal studies post-MI is injectable hyaluronic acid-based hydrogel, which appears to reduce myofiber stress and limit adverse LV remodeling [66, 67].

As a less invasive alternative to intramyocardial injections, biomaterials can also be administered by intracoronary infusion. In 27 post-MI patients, the intracoronary administration of 2 ml of 1% sodium alginate plus 0.3% calcium gluconate within 7 days after MI was feasible and tolerated well [61]. The infusion is assumed to permeate infarcted tissue, where it cross-links into a hydrogel and forms a bioabsorbable cardiac matrix (BCM). The PRESERVATION I randomized controlled trial randomized 303 patients at 2–5 days after ST-elevation MI (STEMI) with TIMI (Thrombolysis In Myocardial Infarction) flow grade 3 in the infarct-related coronary artery to intracoronary application of BCM or saline in a 2:1 ratio [62]. Unfortunately, the investigators did not manage **Fig. 1** An overview of the discussed transcatheter left ventricular restoration devices, grouped by indication for use based on the stage of adverse left ventricular remodeling post myocardial infarction. The represented transvenous annuloplasty device is the mitral loop cerclage, the other devices are differently shaped



to demonstrate a significant difference in adverse clinical event rate or change in LVEDVi from baseline to 6 months, assessed by an independent echocardiography core laboratory blinded to patient data. However, administration of a larger volume of BCM or alternative timing of administration could be worth further investigation. The surgical administration of injectable calcium alginate hydrogel (Algisyl), similar to BCM, has been investigated in patients with advanced heart failure [73]. A total of 78 patients with LVEF $\leq 35\%$, peak VO2 between 9.0 and 14.5 ml/kg/min, and a dilated LV were randomized to treatment with 15 injections or optimal medical therapy. Statistically significant improvements were observed for VO2, 6MWT distance, and NYHA functional class. There was, however, no sham procedure performed in the control group, and there were nine deaths in the Algisyl group (22.5%) vs. four deaths in the control group (10.5%), although the trial was not powered to assess mortality.

Clinical trial design

Study endpoints of clinical trials investigating the effect of LV restoration or enhancing devices should be carefully considered. Because LV dimensions are artificially altered, they should be seen as independent variables. Parameters of LV remodeling such as LVESVi and LVEDVi have to be used

to assess the efficacy of the procedure and should not be interpreted as a proxy variable for patient outcome. Preferable outcome measures to be considered for large pivotal trials are (cardiovascular) mortality, major adverse cardio- and cerebrovascular events, and (non-)heart failure hospitalizations [74]. Other (secondary) endpoints to be considered are functional capacity, severity of dyspnea on visual analogue scale (VAS) or Likert scale, and quality of life assessments. Another important consideration for pivotal randomized controlled trials investigating transcatheter restoration devices is the use of a sham procedure in the control group. A recent meta-analysis suggested that mean improvement in sham groups for percutaneous procedures can be as large as 64%, highlighting the strong placebo effects [75]. A striking example of the importance of a sham procedure is the SYMPLICITY HTN-3 sham controlled trial, which was not able to reproduce the blood pressure-lowering effects of catheter-based renal artery denervation that earlier nonsham controlled trials had observed [76]. A final consideration for device trials is that follow-up after initial results should be extended to 5 years to investigate the long-term effects.

Future perspectives

In May 2017, new European medical device regulations (2017/745) were entered into force [77]. The new regulations increase safety and performance requirements of high-risk devices and the level of supervision of notified bodies. High-risk devices will have to be evaluated by a panel of clinical experts, devices will have to demonstrate equivalence to other safe devices in the market, and more detail will be required in clinical evaluation reports. Relevant information on medical devices will be collected in the European Database on Medical Devices (EUDAMED), and all devices will require a mandatory Unique Device Identification (UDI) so that they can be traced. There will be stricter requirements for postmarket surveillance, including annual safety update reports which have to be assessed by a notified body. There is a transition period until May 2020 during which notified bodies will be designated under the new rules. Before the end of the transition period, all medical devices must be recertified under the new regulations. It will prove to be challenging for small innovative medical device manufacturers to comply with the new requirements.

Patient selection could be very important in determining which patients could have clinical benefit from implantation of a LV restoration device. As the potential benefit of LV restoration at early stages after MI is theoretically larger, early assessment of myocardial viability after MI could be a tool for risk stratification and determining the need for an early intervention. In patients with heart failure, an elementary requirement to be considered for device implantation is that medical therapy has been optimized according to heart failure guidelines, as this has been proven to promote reverse LV remodeling. A potential predictor of poor outcome could be pre-existing diastolic dysfunction, as using a device to artificially reduce the LV volume might further increase LV myocardial stiffness. The lack of consensus on the definition of diastolic dysfunction makes investigating this very challenging. A post hoc analysis of the association between pre-existing diastolic dysfunction and survival in earlier randomized controlled trials on surgical LV restoration could provide more insight.

Improvements in safety and efficacy of mechanical circulatory support devices such as left ventricular assist devices (LVAD) might reduce the need for LV restoration devices in the future. The 6-month results of the HeartMate 3 Left Ventricular Assist System (Abbott, Abbott Park, IL, USA) demonstrate a survival free of any nonsurgical bleeding, thromboembolic event, pump thrombosis, or neurological event in 69% of patients [78]. Progressive improvements could gradually move the indication for use toward earlier stages of heart failure.

Conclusions

Both surgical and transcatheter LV restoration techniques consistently demonstrate improvements in quality-of-life measures and functional status but currently fail to demonstrate a clear survival benefit. Study designs for surgical LV restoration techniques are limited to patients that are already planned for surgery with thoracotomy, such as CABG. The noninvasive nature of transcatheter procedures allows for easier patient selection and identification of independent devicerelated effects. Selection of suitable study endpoints and the use of a sham control procedure is essential in clinical trial design. Transcatheter LV restoration devices show promising results in both animal and in-human studies at different stages of adverse LV remodeling after MI. The Parachute device could be beneficial in heart failure patients with a recent anterior MI, poor systolic function, and a suitable LV anatomy by reducing cardiac dimensions and end-diastolic wall stress. The PARACHUTE IV randomized controlled trial will have to demonstrate additional benefit on top of optimal medical therapy, but it is uncertain whether investigation of the device will be continued. Most of the benefit of LV restoration devices is theoretically gained at early stages after MI, using devices that can alter LV

mechanical properties, such as transcatheter injection of biomaterials in the infarcted region. This will require proper selection of patients at risk of adverse LV remodeling. In the near future, it will prove to be challenging for small medical device manufacturers to comply with requirements of the 2017 European medical device regulations.

Compliance with ethical standards

Conflict of interest The authors declare that they have no conflicts of interest.

Ethical approval This article does not contain any studies with human participants performed by any of the authors.

Open Access This article is distributed under the terms of the Creative Commons Attribution 4.0 International License (http://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided you give appropriate credit to the original author(s) and the source, provide a link to the Creative Commons license, and indicate if changes were made.

References

- Cohn JN, Ferrari R, Sharpe N (2000) Cardiac remodeling—concepts and clinical implications: a consensus paper from an international forum on cardiac remodeling. Behalf of an International Forum on Cardiac Remodeling. J Am Coll Cardiol 35:569–582
- Bhatt AS, Ambrosy AP, Velazquez EJ (2017) Adverse remodeling and reverse remodeling after myocardial infarction. Curr Cardiol Rep 19:71. https://doi.org/10.1007/s11886-017-0876-4
- White HD, Norris RM, Brown MA, Brandt PW, Whitlock RM, Wild CJ (1987) Left ventricular end-systolic volume as the major determinant of survival after recovery from myocardial infarction. Circulation 76:44–51
- 4. Solomon SD, Skali H, Anavekar NS, Bourgoun M, Barvik S, Ghali JK, Warnica JW, Khrakovskaya M, Arnold JM, Schwartz Y, Velazquez EJ, Califf RM, McMurray J, Pfeffer MA (2005) Changes in ventricular size and function in patients treated with valsartan, captopril, or both after myocardial infarction. Circulation 111:3411–3419
- Bolognese L, Neskovic AN, Parodi G, Cerisano G, Buonamici P, Santoro GM, Antoniucci D (2002) Left ventricular remodeling after primary coronary angioplasty: patterns of left ventricular dilation and long-term prognostic implications. Circulation 106:2351–2357
- Moller JE, Hillis GS, Oh JK, Reeder GS, Gersh BJ, Pellikka PA (2006) Wall motion score index and ejection fraction for risk stratification after acute myocardial infarction. Am Heart J 151:419– 425
- Kramer DG, Trikalinos TA, Kent DM, Antonopoulos GV, Konstam MA, Udelson JE (2010) Quantitative evaluation of drug or device effects on ventricular remodeling as predictors of therapeutic effects on mortality in patients with heart failure and reduced ejection fraction: a meta-analytic approach. J Am Coll Cardiol 56:392–406
- Racchi M, Govoni S, Lucchelli A, Capone L, Giovagnoni E (2016) Insights into the definition of terms in European medical device regulation. Expert Rev Med Devices 13:907–917
- Cleutjens JP, Kandala JC, Guarda E, Guntaka RV, Weber KT (1995) Regulation of collagen degradation in the rat myocardium after infarction. J Mol Cell Cardiol 27:1281–1292

- Gajarsa JJ, Kloner RA (2011) Left ventricular remodeling in the post-infarction heart: a review of cellular, molecular mechanisms, and therapeutic modalities. Heart Fail Rev 16:13–21
- Konstam MA, Kramer DG, Patel AR, Maron MS, Udelson JE (2011) Left ventricular remodeling in heart failure: current concepts in clinical significance and assessment. JACC Cardiovasc Imaging 4:98–108
- van den Heuvel AF, van Veldhuisen DJ, van der Wall EE et al (2000) Regional myocardial blood flow reserve impairment and metabolic changes suggesting myocardial ischemia in patients with idiopathic dilated cardiomyopathy. J Am Coll Cardiol 35:19–28
- Yousef ZR, Redwood SR, Marber MS (2000) Postinfarction left ventricular remodelling: where are the theories and trials leading us? Heart 83:76–80
- Springeling T, Kirschbaum SW, Rossi A, Baks T, Karamermer Y, Schulz C, Ouhlous M, Duncker DJ, Moelker A, Krestin GP, Serruys PWJC, de Feyter P, van Geuns RJM (2013) Late cardiac remodeling after primary percutaneous coronary intervention-five-year cardiac magnetic resonance imaging follow-up. Circ J 77:81–88
- Orn S, Manhenke C, Anand IS et al (2007) Effect of left ventricular scar size, location, and transmurality on left ventricular remodeling with healed myocardial infarction. Am J Cardiol 99:1109–1114
- Levine RA, Schwammenthal E (2005) Ischemic mitral regurgitation on the threshold of a solution: from paradoxes to unifying concepts. Circulation 112:745–758
- Writing Group Members, Mozaffarian D, Benjamin EJ et al (2016) Executive summary: heart disease and stroke statistics—2016 update: a report from the American Heart Association. Circulation 133:447–454
- Taniguchi T, Shiomi H, Morimoto T, Watanabe H, Ono K, Shizuta S, Kato T, Saito N, Kaji S, Ando K, Kadota K, Furukawa Y, Nakagawa Y, Horie M, Kimura T (2017) Incidence and prognostic impact of heart failure hospitalization during follow-up after primary percutaneous coronary intervention in ST-segment elevation myocardial infarction. Am J Cardiol 119:1729–1739
- O'Gara PT, Kushner FG, Ascheim DD et al (2013) ACCF/AHA guideline for the management of ST-elevation myocardial infarction: executive summary: a report of the American College of Cardiology Foundation/American Heart Association Task Force on Practice Guidelines. Circulation 127:529–555
- 20. Ibanez B, James S, Agewall S, et al. 2017 ESC guidelines for the management of acute myocardial infarction in patients presenting with ST-segment elevation: the task force for the management of acute myocardial infarction in patients presenting with ST-segment elevation of the European Society of Cardiology (ESC). Eur Heart J 2017
- Kinno M, Nagpal P, Horgan S, Waller AH (2017) Comparison of echocardiography, cardiac magnetic resonance, and computed tomographic imaging for the evaluation of left ventricular myocardial function: part 1 (global assessment). Curr Cardiol Rep 1:19. https:// doi.org/10.1007/s11886-017-0815-4
- 22. Writing Committee Members, Yancy CW, Jessup M et al (2013) ACCF/AHA guideline for the management of heart failure: a report of the American College of Cardiology Foundation/American Heart Association Task Force on Practice Guidelines. Circulation 128:e240–e327
- 23. Yancy CW, Jessup M, Bozkurt B et al (2017) 2017 ACC/AHA/ HFSA focused update of the 2013 ACCF/AHA guideline for the management of heart failure: a report of the American College of Cardiology/American Heart Association Task Force on Clinical Practice Guidelines and the Heart Failure Society of America. Circulation 136:e137–e161
- 24. Ponikowski P, Voors AA, Anker SD et al (2016) ESC guidelines for the diagnosis and treatment of acute and chronic heart failure: the task force for the diagnosis and treatment of acute and chronic heart failure of the European Society of Cardiology (ESC). Developed

with the special contribution of the Heart Failure Association (HFA) of the ESC. Eur J Heart Fail 18:891–975

- Cooley DA, Collins HA, Morris GC Jr, Chapman DW (1958) Ventricular aneurysm after myocardial infarction; surgical excision with use of temporary cardiopulmonary bypass. J Am Med Assoc 167:557–560
- Cooley DA, Frazier OH, Duncan JM, Reul GJ, Krajcer Z (1992) Intracavitary repair of ventricular aneurysm and regional dyskinesia. Ann Surg 215:417–423 discussion 423-4
- Di Donato M, Sabatier M, Dor V, Toso A, Maioli M, Fantini F (1997) Akinetic versus dyskinetic postinfarction scar: relation to surgical outcome in patients undergoing endoventricular circular patch plasty repair. J Am Coll Cardiol 29:1569–1575
- Jones RH, Velazquez EJ, Michler RE, Sopko G, Oh JK, O'Connor CM, Hill JA, Menicanti L, Sadowski Z, Desvigne-Nickens P, Rouleau JL, Lee KL, STICH Hypothesis 2 Investigators (2009) Coronary bypass surgery with or without surgical ventricular reconstruction. N Engl J Med 360:1705–1717
- Carpentier A, Chachques JC (1985) Myocardial substitution with a stimulated skeletal muscle: first successful clinical case. Lancet 1: 1267
- Furnary AP, Jessup FM, Moreira LP (1996) Multicenter trial of dynamic cardiomyoplasty for chronic heart failure. The American Cardiomyoplasty Group. J Am Coll Cardiol 28:1175–1180
- Mann DL, Acker MA, Jessup M, Sabbah HN, Starling RC, Kubo SH, Acorn Trial Principal Investigators and Study Coordinators (2007) Clinical evaluation of the CorCap cardiac support device in patients with dilated cardiomyopathy. Ann Thorac Surg 84: 1226–1235
- 32. Starling RC, Jessup M, Oh JK, Sabbah HN, Acker MA, Mann DL, Kubo SH (2007) Sustained benefits of the CorCap cardiac support device on left ventricular remodeling: three year follow-up results from the Acorn clinical trial. Ann Thorac Surg 84:1236–1242
- 33. Acker MA, Jessup M, Bolling SF, Oh J, Starling RC, Mann DL, Sabbah HN, Shemin R, Kirklin J, Kubo SH (2011) Mitral valve repair in heart failure: five-year follow-up from the mitral valve replacement stratum of the Acorn randomized trial. J Thorac Cardiovasc Surg 142:569–74, 574.e1
- Mann DL, Kubo SH, Sabbah HN, Starling RC, Jessup M, Oh JK, Acker MA (2012) Beneficial effects of the CorCap cardiac support device: five-year results from the Acorn trial. J Thorac Cardiovasc Surg 143:1036–1042
- Wechsler AS, Sadowski J, Kapelak B, Bartus K, Kalinauskas G, Rucinskas K, Samalavicius R, Annest L (2013) Durability of epicardial ventricular restoration without ventriculotomy. Eur J Cardiothorac Surg 44:e189–e192 discussion e192
- Cheng Y, Aboodi MS, Wechsler AS, Kaluza GL, Granada JF, van Bladel K, Annest LS, Yi GH (2013) Epicardial catheter-based ventricular reconstruction: a novel therapy for ischaemic heart failure with anteroapical aneurysm. Interact Cardiovasc Thorac Surg 17: 915–922
- Hernández-Enríquez M, Brugaletta S, Castellá M, et al. New transcatheter treatment of the dilated ischaemic cardiomyopathy with Revivent system. Abstract presented at EuroPCR 2016
- Mazzaferri EL, Jr GS, Sagic D et al (2012) Percutaneous left ventricular partitioning in patients with chronic heart failure and a prior anterior myocardial infarction: results of the PercutAneous Ventricular RestorAtion in Chronic Heart failUre PaTiEnts trial. Am Heart J 163:812–820.e1
- Costa MA, Mazzaferri EL, Jr SH, Abraham WT (2014) Percutaneous ventricular restoration using the parachute device in patients with ischemic heart failure: three-year outcomes of the PARACHUTE first-in-human study. Circ Heart Fail 7:752–758
- Thomas M, Nienaber CA, Ince H, Erglis A, Vukcevic V, Schäfer U, Ferreira RC, Hardt S, Verheye S, Gama Ribeiro V, Sugeng L, Tamburino C (2015) Percutaneous ventricular restoration (PVR)

therapy using the Parachute device in 100 subjects with ischaemic dilated heart failure: one-year primary endpoint results of PARACHUTE III, a European trial. EuroIntervention 11:710–717

- Costa MA, Pencina M, Nikolic S, Engels T, Templin B, Abraham WT (2013) The PARACHUTE IV trial design and rationale: percutaneous ventricular restoration using the parachute device in patients with ischemic heart failure and dilated left ventricles. Am Heart J 165:531–536
- 42. Schmidt T, Frerker C, Thielsen T, Dotz I, Wohlmuth P, Kuck KH, Schäfer U (2014) New evidence for favourable effects on haemodynamics and ventricular performance after Parachute((R)) implantation in humans. Eur J Heart Fail 16:1112–1119
- 43. Yun CH, Sun JY, Templin B, Lin SH, Chen KM, Wu TH, Hung CL, Liu CC, Hsu HH, Alaiti MA, Fares A, DeCicco A, Bezerra HG (2017) Improvements in left ventricular diastolic mechanics after parachute device implantation in patients with ischemia heart failure: a cardiac computerized tomographic study. J Card Fail 23:455– 463
- 44. Tam CC, Fares A, Alaiti A, Shaikh K, Ince H, Erglis A, Bezerra HG, Abraham WT, Costa MA, Attizzani GF (2016) Cardiac computed tomography assessment of the near term impact of percutaneous ventricular restoration therapy (Parachute((R))) on mitral valve geometry. Catheter Cardiovasc Interv 88:E45–E51
- Ladich E, Otsuka F, Virmani R (2014) A pathologic study of explanted parachute devices from seven heart failure patients following percutaneous ventricular restoration. Catheter Cardiovasc Interv 83:619–630
- 46. Grossi EA, Patel N, Woo YJ, Goldberg JD, Schwartz CF, Subramanian V, Feldman T, Bourge R, Baumgartner N, Genco C, Goldman S, Zenati M, Wolfe JA, Mishra YK, Trehan N, Mittal S, Shang S, Mortier TJ, Schweich CJ Jr, RESTOR-MV Study Group (2010) Outcomes of the RESTOR-MV trial (Randomized Evaluation of a Surgical Treatment for Off-Pump Repair of the Mitral Valve). J Am Coll Cardiol 56:1984–1993
- Pedersen WR, Block P, Leon M, Kramer P, Kapadia S, Babaliaros V, Kodali S, Tuzcu EM, Feldman T (2008) iCoapsys mitral valve repair system: percutaneous implantation in an animal model. Catheter Cardiovasc Interv 72:125–131
- Rogers JH, Rahdert DA, Caputo GR, Takeda PA, Palacios IF, Tio FO, Taylor EA, Low RI (2009) Long-term safety and durability of percutaneous septal sinus shortening (the PS(3) system) in an ovine model. Catheter Cardiovasc Interv 73:540–548
- Palacios IF, Condado JA, Brandi S, Rodriguez V, Bosch F, Silva G, Low RI, Rogers JH (2007) Safety and feasibility of acute percutaneous septal sinus shortening: first-in-human experience. Catheter Cardiovasc Interv 69:513–518
- 50. Pighi M, Estevez-Loureiro R, Maisano F, Ussia GP, Dall'Ara G, Franzen O, Laroche C, Settergren M, Winter R, Nickenig G, Gilard M, di Mario C, Transcatheter Valve Treatment Sentinel Registry (TCVT) Investigators of the EURObservational Research Programme (EORP) of the European Society of Cardiology (2017) Immediate and 12-month outcomes of ischemic versus nonischemic functional mitral regurgitation in patients treated with MitraClip (from the 2011 to 2012 pilot sentinel registry of percutaneous edge-to-edge mitral valve repair of the European Society of Cardiology). Am J Cardiol 119:630–637
- 51. Feldman T, Kar S, Elmariah S, Smart SC, Trento A, Siegel RJ, Apruzzese P, Fail P, Rinaldi MJ, Smalling RW, Hermiller JB, Heimansohn D, Gray WA, Grayburn PA, Mack MJ, Lim DS, Ailawadi G, Herrmann HC, Acker MA, Silvestry FE, Foster E, Wang A, Glower DD, Mauri L, EVEREST II Investigators (2015) Randomized comparison of percutaneous repair and surgery for mitral regurgitation: 5-year results of EVEREST II. J Am Coll Cardiol 66:2844–2854
- 52. Harnek J, Webb JG, Kuck KH, Tschope C, Vahanian A, Buller CE, James SK, Tiefenbacher CP, Stone GW (2011) Transcatheter

implantation of the MONARC coronary sinus device for mitral regurgitation: 1-year results from the EVOLUTION phase I study (clinical evaluation of the Edwards Lifesciences percutaneous mitral annuloplasty system for the treatment of mitral regurgitation). JACC Cardiovasc Interv 4:115–122

- 53. Machaalany J, Bilodeau L, Hoffmann R, Sack S, Sievert H, Kautzner J, Hehrlein C, Serruys P, Sénéchal M, Douglas P, Bertrand OF (2013) Treatment of functional mitral valve regurgitation with the permanent percutaneous transvenous mitral annuloplasty system: results of the multicenter international percutaneous transvenous mitral Annuloplasty system to reduce mitral valve regurgitation in patients with heart failure trial. Am Heart J 165:761–769
- Radermecker MA, Lancellotti P, Legrand V, Pierard L (2016) Delayed asymptomatic thrombosis and erosion of the coronary sinus after transvenous mitral annuloplasty with the Viacor device. Ann Thorac Surg 101:2391
- Machaalany J, St-Pierre A, Senechal M et al (2013) Fatal late migration of viacor percutaneous transvenous mitral annuloplasty device resulting in distal coronary venous perforation. Can J Cardiol 29:130.e1–130.e4
- 56. Schofer J, Siminiak T, Haude M, Herrman JP, Vainer J, Wu JC, Levy WC, Mauri L, Feldman T, Kwong RY, Kaye DM, Duffy SJ, Tubler T, Degen H, Brandt MC, van Bibber R, Goldberg S, Reuter DG, Hoppe UC (2009) Percutaneous mitral annuloplasty for functional mitral regurgitation: results of the CARILLON Mitral Annuloplasty Device European Union study. Circulation 120: 326–333
- 57. Siminiak T, Wu JC, Haude M, Hoppe UC, Sadowski J, Lipiecki J, Fajadet J, Shah AM, Feldman T, Kaye DM, Goldberg SL, Levy WC, Solomon SD, Reuter DG (2012) Treatment of functional mitral regurgitation by percutaneous annuloplasty: results of the TITAN trial. Eur J Heart Fail 14:931–938
- Lipiecki J, Siminiak T, Sievert H, et al. Coronary sinus-based percutaneous annuloplasty as treatment for functional mitral regurgitation: the TITAN II trial. Open Heart 2016; 3:e000411-2016-000411. eCollection 2016
- 59. Goldberg SL, Meredith I, Marwick T, Haluska BA, Lipiecki J, Siminiak T, Mehta N, Kaye DM, Sievert H, REDUCE FMR Investigators (2017) A randomized double-blind trial of an interventional device treatment of functional mitral regurgitation in patients with symptomatic congestive heart failure-trial design of the REDUCE FMR study. Am Heart J 188:167–174
- 60. Park YH, Chon MK, Lederman RJ, Sung SC, Je HG, Choo KS, Lee SH, Shin ES, Kim JS, Hwang KW, Lee SY, Chun KJ, Kim CM, Kim JH (2017) Mitral loop cerclage annuloplasty for secondary mitral regurgitation: first human results. JACC Cardiovasc Interv 10:597–610
- Frey N, Linke A, Suselbeck T, Muller-Ehmsen J, Vermeersch P, Schoors D, Rosenberg M, Bea F, Tuvia S, Leor J (2014) Intracoronary delivery of injectable bioabsorbable scaffold (IK-5001) to treat left ventricular remodeling after ST-elevation myocardial infarction: a first-in-man study. Circ Cardiovasc Interv 7: 806–812
- 62. Rao SV, Zeymer U, Douglas PS, al-Khalidi H, Liu J, Gibson CM, Harrison RW, Joseph DS, Heyrman R, Krucoff MW (2015) A randomized, double-blind, placebo-controlled trial to evaluate the safety and effectiveness of intracoronary application of a novel bioabsorbable cardiac matrix for the prevention of ventricular remodeling after large ST-segment elevation myocardial infarction: rationale and design of the PRESERVATION I trial. Am Heart J 170:929–937
- Ryan LP, Matsuzaki K, Noma M, Jackson BM, Eperjesi TJ, Plappert TJ, St. John-Sutton MG, Gorman JH III, Gorman RC (2009) Dermal filler injection: a novel approach for limiting infarct expansion. Ann Thorac Surg 87:148–155

- 64. Dixon JA, Gorman RC, Stroud RE, Mukherjee R, Meyer EC, Baker NL, Morita M, Hamamoto H, Ryan LP, Gorman JH, Spinale FG (2011) Targeted regional injection of biocomposite microspheres alters post-myocardial infarction remodeling and matrix proteolytic pathways. Circulation 124:S35–S45
- 65. McGarvey JR, Kondo N, Witschey WR et al (2015) Injectable microsphere gel progressively improves global ventricular function, regional contractile strain, and mitral regurgitation after myocardial infarction. Ann Thorac Surg 99:597–603
- 66. Rodell CB, Lee ME, Wang H, Takebayashi S, Takayama T, Kawamura T, Arkles JS, Dusaj NN, Dorsey SM, Witschey WRT, Pilla JJ, Gorman JH III, Wenk JF, Burdick JA, Gorman RC (2016) Injectable shear-thinning hydrogels for minimally invasive delivery to infarcted myocardium to limit left ventricular remodeling. Circ Cardiovasc Interv 9:e004058. https://doi.org/10.1161/ CIRCINTERVENTIONS.116.004058
- Dorsey SM, McGarvey JR, Wang H et al (2015) MRI evaluation of injectable hyaluronic acid-based hydrogel therapy to limit ventricular remodeling after myocardial infarction. Biomaterials 69:65–75
- Gooley RP, Meredith IT (2015) The Accucinch transcatheter direct mitral valve annuloplasty system. EuroIntervention 11(Suppl W): W60–W61
- 69. Fieno DS, Kim RJ, Chen EL, Lomasney JW, Klocke FJ, Judd RM (2000) Contrast-enhanced magnetic resonance imaging of myocardium at risk: distinction between reversible and irreversible injury throughout infarct healing. J Am Coll Cardiol 36:1985–1991
- Singelyn JM, Sundaramurthy P, Johnson TD, Schup-Magoffin PJ, Hu DP, Faulk DM, Wang J, Mayle KM, Bartels K, Salvatore M, Kinsey AM, DeMaria AN, Dib N, Christman KL (2012) Catheterdeliverable hydrogel derived from decellularized ventricular extracellular matrix increases endogenous cardiomyocytes and preserves cardiac function post-myocardial infarction. J Am Coll Cardiol 59: 751–763
- 71. Seif-Naraghi SB, Singelyn JM, Salvatore MA, Osborn KG, Wang JJ, Sampat U, Kwan OL, Strachan GM, Wong J, Schup-Magoffin PJ, Braden RL, Bartels K, DeQuach JA, Preul M, Kinsey AM, DeMaria AN, Dib N, Christman KL (2013) Safety and efficacy of an injectable extracellular matrix hydrogel for treating myocardial infarction. Sci Transl Med 5:173ra25
- Wassenaar JW, Gaetani R, Garcia JJ, Braden RL, Luo CG, Huang D, DeMaria AN, Omens JH, Christman KL (2016) Evidence for mechanisms underlying the functional benefits of a myocardial matrix hydrogel for post-MI treatment. J Am Coll Cardiol 67:1074– 1086
- 73. Mann DL, Lee RJ, Coats AJ et al (2016) One-year followup results from AUGMENT-HF: a multicentre randomized controlled clinical trial of the efficacy of left ventricular augmentation with Algisyl in the treatment of heart failure. Eur J Heart Fail 18:314–325
- 74. Zannad F, Garcia AA, Anker SD, Armstrong PW, Calvo G, Cleland JGF, Cohn JN, Dickstein K, Domanski MJ, Ekman I, Filippatos GS, Gheorghiade M, Hernandez AF, Jaarsma T, Koglin J, Konstam M, Kupfer S, Maggioni AP, Mebazaa A, Metra M, Nowack C, Pieske B, Piña IL, Pocock SJ, Ponikowski P, Rosano G, Ruilope LM, Ruschitzka F, Severin T, Solomon S, Stein K, Stockbridge NL, Stough WG, Swedberg K, Tavazzi L, Voors AA, Wasserman SM, Woehrle H, Zalewski A, McMurray JJV (2013) Clinical outcome endpoints in heart failure trials: a European Society of Cardiology Heart Failure Association consensus document. Eur J Heart Fail 15:1082–1094
- 75. Jonas WB, Crawford C, Colloca L, et al. To what extent are surgery and invasive procedures effective beyond a placebo response? A systematic review with meta-analysis of randomised, sham controlled trials. BMJ Open 2015; 5:e009655-2015-009655
- 76. Bhatt DL, Kandzari DE, O'Neill WW, D'Agostino R, Flack JM, Katzen BT, Leon MB, Liu M, Mauri L, Negoita M, Cohen SA,

Oparil S, Rocha-Singh K, Townsend RR, Bakris GL, SYMPLICITY HTN-3 Investigators (2014) A controlled trial of renal denervation for resistant hypertension. N Engl J Med 370: 1393–1401

- Official Journal of the European Union. Regulation (EU) 2017/745 of the European Parliament and of the Council of 5 April 2017 on medical devices, amending Directive 2001/83/EC, Regulation (EC) No 178/2002 and Regulation (EC) No 1223/2009 and repealing Council Directives 90/385/EEC and 93/42/EEC. 2017:L 117
- 78. Uriel N, Colombo PC, Cleveland JC, Long JW, Salerno C, Goldstein DJ, Patel CB, Ewald GA, Tatooles AJ, Silvestry SC, John R, Caldeira C, Jeevanandam V, Boyle AJ, Sundareswaran KS, Sood P, Mehra MR (2017) Hemocompatibility-related outcomes in the MOMENTUM 3 trial at 6 months: a randomized controlled study of a fully magnetically levitated pump in advanced heart failure. Circulation 135:2003–2012